

The Healing Touch



*It begins with products
and services from
National Wound Care.*



Because your patients count on you, you can count on National Wound Care.

As a health care professional, you gently care for the needs of each wound care patient.

As your partner, National Wound Care provides the products and services you need to deliver that healing touch.

And our ability to do just that has established National Wound Care as a leader in providing the highest quality wound care equipment. Nursing homes, home health agencies, hospice agencies, and hospitals across the country turn to us to assist them in treating and preventing wounds. They know they can count on us to deliver the most innovative, cost-effective solutions for their wound care patients.

We make sure you have state-of-the-art tools to assist you—from the latest in support surface to mattress replacements to wheelchair cushions. And when you need more than tools, the highly trained professionals at National Wound Care are just a phone call away. Whether you have questions about products or support services or comprehensive wound care management, we can help.

When you choose National Wound Care, you get superior products and much more. You get the peace of mind that comes with knowing you have a partner in delivering the finest care. We invite you to review our complete line of wound care products and programs.



**National Wound Care
Closed Suction Wound
Drainage System**
Code name: *Wound Sucker™*

The National Wound Care Closed Suction Wound Drainage System (CSWDS) is a powered suction machine connected via tubing directly into the patient's wound bed used to remove exudate from the wound and to apply a negative pressure to the wound site. The NWC system utilizes a special suction pump which has FDA approval for use in wound care for wound drainage, and because of the effect of having an 'airtight seal' over the wound, in a 'closed suction wound drainage' application, it results in a 'negative pressure' applied to the wound site. Illinois Department of Public Aid has approved the NWC CSWDS under the "Negative Pressure Wound Therapy Pump", procedure code E2402. For more information on this new product contact your local NWC representative, or contact our office at 1-800-982-1835.



Here to Serve

It is not just the superior products that set NWC above our competition, it is our dedicated and compassionate staff that are the core of who we are. Every individual that bears the NWC name has been hired because of possessing a servant's heart, and a commitment and desire to serve our customers and clients with the kindness and compassion by which they want to be served. Take time to read through our *Core Values* on the back of this brochure, as they were developed by all the individuals that comprise *National Wound Care*.



NWC 5000 Kinetic Low Air-Loss

A dynamic, computer-controlled low air loss mattress replacement system with kinetic features. Converts any hospital bed into a pressure-relief surface. Pressure relief coupled with two dynamic modalities (derma-wave and derma-pulse) increases capillary blood flow, which in turn accelerates healing.



NWC 4500 Low Air-Loss w/Lateral Rotation

Combines kinetic rotation with the convenience of a portable low air loss mattress replacement system. Turns a patient 40 degrees bilaterally, reducing respiratory complications associated with immobility.




NWC 4000 Low Air-Loss

A low air loss mattress replacement system that converts any hospital bed into an effective pressure reduction therapy surface. The 9-inch depth affords optimum weight distribution as well as enhanced patient comfort.



NWC 3100 Alternating Pressure Mattress

A mattress replacement system that provides superior pressure relief through the use of deep cell alternating pressure therapy.




NWC 2000ET Self Adjusting Mattress

Non-powered, self adjusting pressure relief support system.



NWC 1000 Alternating Pressure Overlay

Combines the portability of an overlay with the therapeutic effectiveness of alternating pressure relief.



NWC 400 Gel Overlay

Our most effective Group 1 support surface, this lightweight, static overlay is constructed of a solid gel core surrounded by a foam layer for superior comfort and pressure reduction.



Preventive Products Foam Mattresses and Wheelchair Cushions

National Wound Care offers additional products to help prevent skin breakdown for at-risk patients. Contact us for more information on our complete line of pressure reduction mattress replacement systems and our special lightweight wheelchair cushions.

Bariatric Products

Bed/mattress systems also available for bariatric patients up to 1000 lbs.



Core Values

- We strive to deal honestly and openly with each other and our customers, and conduct our business with integrity.
- We believe in demonstrating kindness and compassion in all our dealings.
- We will treat others with the courtesy and respect by which we wish to be treated.
- We strive to excel in all areas of responsibility.
- We will display teamwork in all endeavors.
- We strive to enhance communication both internally and externally.
- We will be adaptable to the changing needs of our business.
- We will encourage and promote the growth and development of each employee.
- We will encourage each individual to be the best he/she can be.
- We are here to serve our customers and one another.

Mission Statement

Our mission is to be the preferred company of each customer we serve. Our principle focus is to enhance the quality of life of the patient by providing our clients with comprehensive wound care equipment that meet or exceed their expectations. We are committed to the success of our customers by assisting caregivers in treating and preventing wounds in the most effective, cost-efficient manner.



2906 W. Clark Drive
Champaign, Illinois 61822
1-800-982-1835